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ESPAÑOL

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INSTRUCTIONS MANUAL FOR SERVICE MAINTENANCE USE

www.kmdargentina.com.ar

CONGRATULATIONS!

You have joined the KMD technology user group. To ensure a long use and optimal performance of your equipment, please read this instructions manual carefully.

FEATURES

- > Light-alloy metallic body
- > Anti-glare chromium treatment
- > Autoclavable up to 135°C
- > Spray throughout the work area
- > Precision metallic axle

PLANT AND ADMINISTRATION

José Ignacio Rucci 3935/39 (C1439FVQ)
CABA - Buenos Aires - Argentina
Tel./Fax: 54-11 4601-3208 / 4638-8282

Manufactured in Argentina
by Meca Dent S.R.L.


| HIGH-SPEED MODELS | SPEED RPM | COUPLING | MAX. AIR PRESSURE | SPRAY | WATER PRESSURE | CONSUMPTION AIR | CONSUMPTION WATER | FASTENING CLAMP | WEIGHT Gms. | BURR TYPE FG Standard 01.6 -0.01 mm | AUTOCLAVE STERILIZATION 135°C(275°F) H2O Steam |
|---------------------------|-------------|----------|--------------------|----------|-------------------|-----------------|-------------------|-----------------|-------------|-------------------------------------|--|
| Maxtorque | 300000 ±10% | X X | 2,1 kg/cm² (30PSI) | Simple | 250 Kpa (2.5 bar) | Max. 66 N/min | Min. 50ml/min | | | | |
| Standard | 270000 ±10% | X | 2,8 kg/cm² (40PSI) | Quaduple | X | X | X | FG Chuck | 45 | X | X |
| S2000 | 300000 ±10% | X | | X | X | X | X | FG Chuck | 45 | X | X |
| Belair LED* | 375000 ±8% | | X X | X | X | X | X | PB Chuck | 45 | X | X |
| Belair | 280000 ±10% | X X | X | X | X | X | X | PB Chuck | 45 | X | X |
| Brio (Ultra torque) | 380000 ±8% | X | X | X | X | X | X | PB Chuck | 45 | X | X |
| Colibri NS | 330000 ±10% | X X | X | X | X | X | X | PB Chuck | 45 | X | X |
| S2000 | 300000 ±10% | X X | X | X | X | X | X | PB Chuck | 45 | X | X |
| SIRION -AGR-Multiflex* | 295000 ±5% | | X | X | X | X | X | PB Chuck | 65 | X | X |
| Colibri NS -AGR-Multiflex | 350000 ±6% | | X | X | X | X | X | PB Chuck | 60 | X | X |
| LABns | 300000 ±10% | X X | X | X | X | X | X | FG Chuck | 70 | X | X |

*Built-in lighting means

INSTALLATION AND OPERATION

Check whether the equipment's coupling is compatible with the handpiece's coupling, and then check whether air and water are coming out through the connector, screw the coupling nut clockwise and tighten. (*For Multiflex coupling, connect the turbine to the coupling until a "click" is heard.*) Upon completing this operation, the installation is finished. *It is essential to provide clean and dry compressed air.* For the handpiece's correct functioning, a maximum air pressure of (see table) is required. Water flow regulation: With the equipment switch in the open position and the handpiece in operation, close it gradually until obtaining the desired spray. To place and change the burrs correctly, see pictures (according to the model). The burr must always reach the end to be fixed correctly. *Use high quality and standard burrs.*

MODELS WITH FG CHUCK

DISINFECTION AND STERILIZATION

Decouple the equipment handpiece. Clean the outer part with alcohol. Do not immerse in disinfectant solutions. Inject KMD spray oil for 2 seconds. See pictures (according to coupling). Wrap in paper. Sterilize in autoclave up to 135°C (275°F) during one cycle. Remove the handpiece and couple it with the air conditioning line, running it for a few seconds for the excess oil to be eliminated.

STATEMENT: Both the burr remover and the burr key for the Maxtorque, Standard, S2000 and LABns models are sterilizable by autoclave.

LUBRICATION

2-WAY, 4-WAY AND H6 COUPLING

HOW TO USE

MULTIFLEX COUPLING

MAINTENANCE

Lubrication is required after each work day. We recommend the use of KMD spray oil by injecting it for 2 seconds, following the above pictures. In the event of very hard work, lubrication is required during the work day and at the end of the day.

IMPORTANT: IF YOU HAVE ANY TROUBLE, PLEASE CHECK THIS TABLE BEFORE SENDING THE EQUIPMENT FOR TECHNICAL SERVICE:

FAILURE

CAUSE

PRECAUTION

NEVER immerse the handpiece in disinfectant or cleaning solutions. DO NOT sterilize before removing the burr. In PB models, do not press the button while the piece is in motion.

The piece's number of revolutions decreases suddenly (it is acoustically perceived as a change in tone)

The spray refrigeration does not work.

The burr goes out during operation.

The burr does not fit.

1 - There is water in the air line. The anchor between the piece and the hose is loose. 2 - Defective piece hose. 3 - Decrease in air flow. 4 - The rotor is broken down. The handpiece lid is loose.

The water supply is shut off.

Worn-out burr shaft or poor quality burr.

Burr shaft out of tolerance (Ø max. 1.60 mm).

1 - Fasten the coupling securely with the piece in operation. 2 - Change the hose. 3 - Check that air is supplied to the coupling. 4 - Send the equipment for technical service.

Check that water is supplied from the equipment.

Carefully clear the spray holes with a probe and inject pressure air from the head to the equipment's coupling.

Replace with a new high-quality burr.

Replace with a new high-quality burr.

REV. 1 -02/15

INITIAL MAINTENANCE

Clean, inject oil in, and sterilize the handpiece before using it for the first time. Follow the disinfection and sterilization procedure.